

Anti-Plzf Polyclonal Antibody For Wb, Ihc-P, Elisa - Q05516

Item Number: CM0030271



Introduction

Rabbit polyclonal antibody against PLZF (ZBTB16/ZNF145) validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Suitable for studying transcriptional repression and myeloid differentiation.

[Learn More](#)

Parameter	Value
Product Name	PLZF Rabbit pAb
Remarks/Alias	PLZF; ZNF145
Species	Human
Gene ID	7704
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 199-300 of human PLZF (NP_005997.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05516
Protein Weight	74kDa
Shipping	Ice bag