

Anti-Plzf Rabbit Monoclonal Antibody For Wb - Q05516

Item Number: CM0002570



Introduction

High-specificity rabbit monoclonal antibody targeting human PLZF (ZBTB16/ZNF145). Validated for WB, IP, IF-P, ELISA. Cross-reacts with mouse and rat. Ideal for transcriptional regulation and ubiquitination studies.

[Learn More](#)

Parameter	Value
Product Name	PLZF Rabbit mAb
Remarks/Alias	PLZF; ZNF145
Species	Human
Gene ID	7704
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 199-300 of human PLZF (NP_005997.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05516
Protein Weight	74kDa
Shipping	Ice bag