


Anti-Pmaip1 (Noxa/Apr) Polyclonal Antibody For Wb, If/Icc, Elisa - Q13794

Item Number: CM0025355



Introduction

Rabbit polyclonal antibody against human PMAIP1 (NOXA), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for apoptosis and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	PMAIP1 Rabbit pAb
Remarks/Alias	APR, NOXA, PMAIP1
Species	Human
GeneID (Human)	5366
GeneID	5366
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-54 of human PMAIP1 (NP_066950.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13794
Protein Weight	6kDa
Shipping	Ice bag