

Anti-Pmch Rabbit Polyclonal Antibody For Wb, Elisa - P20382

Item Number: CM0026030



Introduction

Rabbit polyclonal antibody against human PMCH, validated for Western blot (WB) and ELISA. Reacts with human and mouse samples, targeting the pro-MCH precursor protein (~19 kDa). Ideal for neuropeptide research, obesity, and metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	PMCH Rabbit pAb
Remarks/Alias	MCH; ppMCH; PMCH
Species	Human
Gene ID (Human)	5367
Gene ID	5367
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-165 of human PMCH (NP_002665.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P20382
Protein Weight	19kDa
Shipping	Ice bag