

Anti-Pmel (Gp100) Monoclonal Antibody For Wb, Ihc-P, Elisa - P40967

Item Number: CM0004872

Introduction

Rabbit monoclonal antibody against human PMEL (gp100), validated for Western blot, IHC (paraffin), and ELISA. Suitable for melanoma research, melanosome biology, and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	PMEL Rabbit mAb
Remarks/Alias	P1; SI; SIL; ME20; P100; SILV; HMB45; ME20M; gp100; HMB-45; ME20-M; PMEL17; D12S53E; PMEL
Species	Human
GeneID	6490
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P40967
Protein Weight	70kDa
Shipping	Ice bag