

Anti-Pmel Polyclonal Antibody For Wb, If/Icc, Elisa - P40967

Item Number: CM0027855



Introduction

Rabbit polyclonal antibody targeting human PMEL (gp100), validated for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat, ideal for melanosome biogenesis and melanoma antigen studies. 70 kDa.

[Learn More](#)

Parameter	Value
Product Name	PMEL Rabbit pAb
Remarks/Alias	P1; SI; SIL; ME20; P100; SILV; HMB45; ME20M; gp100; HMB-45; ME20-M; PMEL17; D12S53E; PMEL
Species	Human
GeneID (Human)	6490
GeneID	6490
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 320-430 of human PMEL (NP_008859.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40967
Protein Weight	70kDa
Shipping	Ice bag