

Anti-Pmepa1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q969W9

Item Number: CM0021429

Introduction

High-quality rabbit polyclonal antibody targeting human PMEPA1 (TMEPAI) protein, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for TGF-beta signaling and prostate cancer research.

[Learn More](#)

Parameter	Value
Product Name	PMEPA1 Rabbit pAb
Remarks/Alias	STAG1; TMEPAI; PMEPA1
Species	Human
GenelD (human)	56937
GenelD	56937
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 173-252 of human PMEPA1 (NP_954638.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969W9
Protein Weight	32kDa
Shipping	Ice bag