

Anti-Pmf1 Monoclonal Antibody For Wb, Elisa - Q6P1K2

Item Number: CM0003184



Introduction

Rabbit monoclonal antibody against PMF1, suitable for Western blot and ELISA. Validated in human and rat samples. 23kDa polyamine-modulated factor 1 involved in chromosome segregation and gene transcription.

[Learn More](#)

Parameter	Value
Product Name	PMF1 Rabbit mAb
Remarks/Alias	PMF1
Species	Human
Gene ID (Human)	11243
Gene ID	11243
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human PMF1 (NP_009152.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q6P1K2
Protein Weight	23kDa
Shipping	Ice bag