

Anti-Pml Polyclonal Antibody For Wb, Ihc-P, Elisa - P29590

Item Number: CM0024152



Introduction

Rabbit polyclonal antibody against human PML protein, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. PML is a scaffold protein of nuclear bodies involved in tumor suppression, apoptosis, and antiviral defense. Ideal for IVD research use.

[Learn More](#)

Parameter	Value
Product Name	PML Rabbit pAb
Remarks/Alias	MYL; RNF71; PP8675; TRIM19; PML
Species	Human
GeneID (Human)	5371
GeneID	5371
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 571-611 of human PML (NP_150241.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29590
Protein Weight	98kDa
Shipping	Ice bag