

Anti-Pml+Rara Rabbit Polyclonal Antibody For If/Icc And Elisa - P29590 / P10276

Item Number: CM0020099



Introduction

Polyclonal antibody targeting the PML-RARA fusion protein (type L APL). Validated for IF/ICC and ELISA in human and rat samples. Ideal for APL leukemia research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	PML+RARA Rabbit pAb
Remarks/Alias	PML/RARA; Type L APL
Species	Human
GeneID	5371/5914
Immunogen	synthetic peptide corresponding to amino acids 500-600/1-100 of human PML+RARA (NP_150241.2/NP_000955.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P29590/P10276
Protein Weight	47-48kDa/62-97kDa/39kDa/50kDa
Shipping	Ice bag