

Anti-Pmm2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O15305

Item Number: CM0025400



Introduction

Rabbit polyclonal antibody against human PMM2 (Phosphomannomutase 2), suitable for Western blot, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat samples. Useful for congenital disorders of glycosylation research.

[Learn More](#)

Parameter	Value
Product Name	PMM2 Rabbit pAb
Remarks/Alias	PMI; CDG1; CDGS; PMI1; CDG1a; PMM 2; PMM2
Species	Human
GeneID (Human)	5373
GeneID	5373
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-246 of human PMM2 (NP_000294.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15305
Protein Weight	28kDa
Shipping	Ice bag