

Anti-Pmp70/Abcd3 Monoclonal Antibody For Wb, If/Icc, Elisa - P28288

Item Number: CM0011901



Introduction

Rabbit monoclonal antibody against human PMP70/ABCD3 (P28288), validated for WB, IF/ICC, ELISA with cross-reactivity to mouse/rat. Ideal for peroxisomal fatty acid transport and beta-oxidation research.

[Learn More](#)

Parameter	Value
Product Name	PMP70/ABCD3 Rabbit mAb
Remarks/Alias	ZWS2; ABC43; CBA55; PMP70; PXMP1; PMP70/ABCD3
Species	Human
Gene ID (Human)	5825
Gene ID	5825
Immunogen	Synthetic peptide corresponding to amino acids 560-659 of human PMP70/ABCD3 (P28288)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28288
Protein Weight	75kDa
Shipping	Ice bag