

Anti-Pmp70/Abcd3 Rabbit Polyclonal Antibody - P28288

Item Number: CM0023376



Introduction

High-quality rabbit polyclonal antibody against PMP70/ABCD3, validated for WB, IF/ICC, IHC-P, ELISA. Suitable for studying peroxisomal fatty acid transport and metabolism in human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	PMP70/ABCD3 Rabbit pAb
Remarks/Alias	ZWS2; ABC43; CBAS5; PMP70; PXMP1; PMP70/ABCD3
Species	Human
GeneID (human)	5825
GeneID	5825
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-475 of human PMP70/ABCD3 (NP_002849.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28288
Protein Weight	75 kDa
Shipping	Ice bag