

Anti-Pmpca Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Ip, Elisa - Q10713

Item Number: CM0006729

Introduction

High-specificity rabbit monoclonal antibody targeting PMPCA (alpha-MPP). Validated for WB, IF/ICC, IHC-P, IP, ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial processing studies.

[Learn More](#)

Parameter	Value
Product Name	PMPCA Rabbit mAb
Remarks/Alias	CLA1; CPD3; MAS2; P-55; SCAR2; INPP5E; Alpha-MPP; PMPCA
Species	Human
GeneID (human)	23203
GeneID	23203
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 291-525 of human PMPCA (NP_055975.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q10713
Protein Weight	58 kDa
Shipping	Ice bag