

Anti-Pmpcb Rabbit Polyclonal Antibody For Wb, Elisa - O75439

Item Number: CM0015635

Introduction

High-quality rabbit polyclonal antibody targeting PMPCB, a key mitochondrial processing protease. Validated for WB and ELISA applications in human, mouse, and rat samples. Ideal for mitochondrial import and PINK1 turnover research.

[Learn More](#)

Parameter	Details
Product Name	PMPCB Rabbit pAb
Remarks/Alias	MAS1; MPPB; P-52; MPP11; MPPP52; Beta-MPP; PMPCB
Species	Human
GeneID (Human)	9512
GeneID	9512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-380 of human PMPCB (NP_004270.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75439
Protein Weight	54kDa
Shipping	Ice bag