

Anti-Pms1 Rabbit Polyclonal Antibody For Wb And Elisa - P54277

Item Number: CM0015842



Introduction

Rabbit polyclonal PMS1 antibody for WB and ELISA in human, mouse, rat. Targets DNA mismatch repair protein PMS1, linked to hereditary nonpolyposis colorectal cancer.

[Learn More](#)

Parameter	Value
Product Name	PMS1 Rabbit pAb
Remarks/Alias	MLH2; PMSL1; hPMS1; HNPCC3; PMS1
Species	Human
GeneID (Human)	5378
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 405-619 of human PMS1 (NP_000525.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P54277
Protein Weight	106kDa
Shipping	Ice bag