


Anti-Pms1 Rabbit Polyclonal Antibody For Wb, Elisa - P54277

Item Number: CM0021482



Introduction

High-quality rabbit polyclonal antibody targeting human PMS1 (DNA mismatch repair protein). Validated for WB and ELISA applications, supporting research in DNA repair and hereditary nonpolyposis colorectal cancer (HNPCC).

[Learn More](#)

Parameter	Value
Product Name	PMS1 Rabbit pAb
Remarks/Alias	MLH2; PMSL1; hPMS1; HNPCC3; PMS1
Species	Human
GeneID (Human)	5378
GeneID	5378
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 405-619 of human PMS1 (NP_000525.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P54277
Protein Weight	106kDa
Shipping	Ice bag