


Anti-Pms2 Rabbit Monoclonal Antibody (Ko Validated) For Wb And Elisa - P54278

Item Number: CM0006435



Introduction

Anti-PMS2 Rabbit Monoclonal Antibody (KO Validated) optimized for Western blot and ELISA, recognizing human PMS2 (96 kDa). A mismatch repair endonuclease marker for colorectal cancer (HNPCC) research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PMS2 Rabbit mAb
Remarks/Alias	MLH4; PMS-2; PMSL2; HNPCC4; LYNCH4; MMRCS4; PMS2CL; S2
Species	Human
GeneID (Human)	5395
Immunogen	Recombinant protein of human PMS2
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P54278
Protein Weight	96kDa
Shipping	Ice bag