


Anti-Pms2 Rabbit Polyclonal Antibody For Wb, Elisa - P54278

Item Number: CM0025083



Introduction

Rabbit polyclonal antibody targeting human PMS2, a mismatch repair endonuclease. Validated for WB and ELISA; cross-reacts with rat. Ideal for cancer research and MMR studies.

[Learn More](#)

Parameter	Value
Product Name	PMS2 Rabbit pAb
Remarks/Alias	MLH4; PMS-2; PMSL2; HNPCC4; LYNCH4; MMRC54; PMS2CL; PMS2
Species	Human
GeneID (Human)	5395
GeneID	5395
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P54278
Protein Weight	96kDa
Shipping	Ice bag