


Anti-Pms2 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P54278

Item Number: CM0013463



Introduction

Rabbit polyclonal antibody targeting human PMS2 mismatch repair endonuclease, validated for WB, IF/ICC, IP, ELISA. Cross-reacts with mouse PMS2. Ideal for cancer and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PMS2 Rabbit pAb
Remarks/Alias	MLH4; PMS-2; PMSL2; HNPCC4; LYNCH4; MMRC4; PMS2CL; S2
Species	Human
GeneID	5395
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 390-670 of human PMS2 (NP_000526.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P54278
Protein Weight	96kDa
Shipping	Ice bag