


Anti-Pmvk Polyclonal Antibody For Wb - Q15126

Item Number: CM0013370



Introduction

Rabbit polyclonal antibody against human PMVK (Phosphomevalonate kinase), validated for WB and ELISA, cross-reacts with mouse and rat. Key reagent for mevalonate pathway and cholesterol biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	PMVK Rabbit pAb
Remarks/Alias	PMK; PMKA; PMKASE; POROK1; HUMPMKI; PMVK
Species	Human
GeneID (Human)	10654
GeneID	10654
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-192 of human PMVK (NP_006547.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15126
Protein Weight	22kDa
Shipping	Ice bag