


Anti-Pnkp Rabbit Monoclonal Antibody For Wb, Elisa - Q96T60

Item Number: CM0006549



Introduction

Anti-PNKP rabbit monoclonal antibody validated for Western blotting (WB) and ELISA. Cross-reacts with human, mouse, and rat. Targets the DNA repair protein PNKP, ideal for BER and NHEJ pathway studies.

[Learn More](#)

Parameter	Value
Product Name	PNKP Rabbit mAb
Remarks/Alias	PNK; AOA4; MCSZ; CMT2B2; EIEE10; PNKP
Species	Human
GeneID (Human)	11284
GeneID	11284
Immunogen	Synthetic peptide (A synthetic peptide corresponding to a sequence within amino acids 400-500 of Human PNKP (NP_009185.2).)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96T60
Protein Weight	57kDa
Shipping	Ice bag