


Anti-Pnkp Rabbit Polyclonal Antibody For Wb - Q96T60

Item Number: CM0028500



Introduction

PNKP Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat PNKP, a bifunctional kinase/phosphatase essential for DNA repair. Ideal for DNA damage repair studies.

[Learn More](#)

Parameter	Value
Product Name	PNKP Rabbit pAb
Remarks/Alias	PNK; AOA4; MCSZ; CMT2B2; EIEE10; PNKP
Species	Human
Gene ID (Human)	11284
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-320 of human PNKP (NP_009185.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96T60
Protein Weight	57kDa
Shipping	Ice bag