

## Anti-Pnliprp1 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P54315

Item Number: CM0020210



### Introduction

Rabbit polyclonal antibody targeting human inactive pancreatic lipase-related protein 1 (PNLIPRP1/PLRP1), validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat, ideal for pancreatic research.

[Learn More](#)

| Parameter        | Specification                                                                                                         |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | PNLIPRP1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | PLRP1; PNLIPRP1                                                                                                       |
| Species          | Human                                                                                                                 |
| GeneID (Human)   | 5407                                                                                                                  |
| GeneID           | 5407                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 298-467 of human PNLIPRP1 (NP_006220.1) |
| Source           | Rabbit                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                 |
| Application      | WB, IF-P, IHC-P, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                     |
| SWISS            | P54315                                                                                                                |
| Protein Weight   | 52kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |