

Anti-Pnliprp2 Polyclonal Antibody For Wb, If-P, Elisa - P54317

Item Number: CM0015553



Introduction

High-quality rabbit polyclonal antibody targeting human PNLIPRP2/PLRP2. Validated for WB, IF-P, ELISA. Cross-reacts with mouse and rat. Ideal for lipid metabolism, pancreatic digestion, and neuronal studies.

[Learn More](#)

Parameter	Value
Product Name	PNLIPRP2 Rabbit pAb
Remarks/Alias	PLRP2; PNLIPRP2
Species	Human
Gene ID (Human)	5408
Gene ID	5408
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 301-470 of human PNLIPRP2 (NP_005387.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P54317
Protein Weight	52kDa
Shipping	Ice bag