


Anti-Pnma2 Polyclonal Antibody For Wb And Elisa - Q9UI42

Item Number: CM0031047



Introduction

Rabbit polyclonal antibody targeting human PNMA2 (paraneoplastic Ma2 antigen), validated for Western blotting and ELISA. Suitable for neuro-oncology and paraneoplastic syndrome research.

[Learn More](#)

Parameter	Value
Product Name	PNMA2 Rabbit pAb
Alias/Remarks	MA2; MM2; RGAG2; PNMA2
Species	Human
Gene ID	10687
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human PNMA2 (NP_009188.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS Accession	Q9UL42
Protein Weight	42kDa
Shipping	Ice bag