


Anti-Pno1 Polyclonal Antibody For Wb, Elisa - Q9Nrx1

Item Number: CM0022410



Introduction

Rabbit polyclonal antibody targeting human PNO1, a nucleolar RNA-binding protein involved in ribosome biogenesis. Validated for Western blot and ELISA applications with cross-reactivity to human and rat samples. Ideal for ribosomal research.

[Learn More](#)

Parameter	Value
Product Name	PNO1 Rabbit pAb
Remarks/Alias	RRP20; KHRBP1; PNO1
Species	Human
GeneID (Human)	56902
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-252 of human PNO1 (NP_064528.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9NRX1
Protein Weight	28kDa
Shipping	Ice bag