

Anti-Pnpo Rabbit Polyclonal Antibody For Wb - Q9Nvs9

Item Number: CM0020035



Introduction

This PNPO antibody is validated for Western blot, immunofluorescence, and ELISA. It reacts with human, mouse, and rat samples. Use for research on vitamin B6 metabolism.

[Learn More](#)

Parameter	Specification
Product Name	PNPO Rabbit pAb
Remarks/Alias	PDXPO; HEL-S-302; PNPO
Species	Human
GeneID (Human)	55163
GeneID	55163
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-261 of human PNPO (NP_060599.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NVS9
Protein Weight	30kDa
Shipping	Ice bag