

Anti-Poc1A Polyclonal Antibody For Wb, Elisa - Q8Nbt0

Item Number: CM0029499

Introduction

Anti-POC1A rabbit polyclonal antibody for WB and ELISA; cross-reacts with human, mouse, rat. Targets centriolar protein POC1A involved in centriole assembly and ciliogenesis.

[Learn More](#)

Parameter	Value
Product Name	POC1A Rabbit pAb
Remarks/Alias	PIX2; SOFT; WDR51A; POC1A
Species	Human
GeneID (Human)	25886
GeneID	25886
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 138-407 of human POC1A (NP_056241.3))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NBT0
Protein Weight	45kDa
Shipping	Ice bag