


Anti-Podoplanin Rabbit Monoclonal Antibody For Wb, Elisa - Q86YI7

Item Number: CM0005206



Introduction

Anti-Podoplanin rabbit monoclonal antibody (UniProt Q86YL7) for WB and ELISA. Detects human, mouse, and rat Podoplanin. Useful in cancer, cell migration, and adhesion studies.

[Learn More](#)

Parameter	Specification
Product Name	Podoplanin Rabbit mAb
Remarks/Alias	T1A; GP36; GP40; Gp38; OTS8; T1A2; T1A; D2-40; T1A-2; AGGRUS; HT1A-1; PA2.26; Podoplanin
Species	Human
GeneID (human)	10630
GeneID	10630
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 63-162 of human Podoplanin (Q86YL7).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86YL7
Protein Weight	17 kDa
Shipping	Ice bag