

Anti-Podoplanin Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q62011

Item Number: CM0011384



Introduction

Validated for WB, IF-P, IHC-P, and ELISA, this rabbit monoclonal antibody targets mouse Podoplanin (Gp38/T1alpha). Ideal for studies on lymphatic development, cancer metastasis, and cell migration. Specific to mouse, ~18 kDa.

[Learn More](#)

Parameter	Value
Product Name	Podoplanin Rabbit mAb
Remarks/Alias	E11; T1a; Gp38; OTS-8; T1alpha; RANDAM-2; T1-alpha
Species	Mouse
GeneID (Human)	10630
GeneID	14726
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-133 of mouse Podoplanin (NP_034359.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	Q62011
Protein Weight	18 kDa
Shipping	Ice bag