

Anti-Podoplanin Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q86YI7

Item Number: CM0013473



Introduction

Rabbit polyclonal antibody against human Podoplanin (PDPN, Aggrus, GP36), validated in WB, IHC-P, ELISA; reacts with human, mouse, rat. For tumor migration, lymphatic and platelet research.

[Learn More](#)

Parameter	Value
Product Name	Podoplanin Rabbit pAb
Remarks/Alias	T1A; GP36; GP40; Gp38; OTS8; T1A2; T11A; D2-40; T1A-2; AGGRUS; HT1A-1; PA2.26; Podoplanin
Species	Human
Gene ID (Human)	10630
Gene ID	10630
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-131 of human Podoplanin (NP_006465.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86YL7
Protein Weight	17kDa
Shipping	Ice bag