

Anti-Podxl Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O00592

Item Number: CM0009218

Introduction

Rabbit monoclonal antibody against human PODXL, validated for Western blot, immunohistochemistry (paraffin), and ELISA. Suitable for studying cell adhesion, podocyte biology, and cancer progression. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	PODXL Rabbit mAb
Remarks/Alias	PC; PDX; PCLP; Gp200; gp135; PCLP-1; PODXL1; PODXL
Species	Human
GeneID (Human)	5420
GeneID	5420
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-479 of human PODXL (O00592).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O00592
Protein Weight	59kDa
Shipping	Ice bag