


Anti-Pof1B Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wvv4

Item Number: CM0023251



Introduction

POF1B Rabbit Polyclonal Antibody validated for WB and ELISA, reactive to human and rat. Suitable for research on epithelial monolayer organization and ovarian development. SWISS: Q8WVV4.

[Learn More](#)

Parameter	Value
Product Name	POF1B Rabbit pAb
Remarks/Alias	POF; POF2B; POF1B
Species	Human
GeneID (Human)	79983
GeneID	79983
Immunogen	Recombinant protein, Recombinant fusion protein containing a sequence corresponding to amino acids 330-560 of human POF1B (NP_079197.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8WVV4
Protein Weight	68kDa
Shipping	Ice bag