

Anti-Poglut1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nbl1

Item Number: CM0022848



Introduction

Rabbit polyclonal antibody targeting human POGLUT1 (Rumi, CLP46), validated in WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for Notch signaling, glycosylation, and muscle development studies.

[Learn More](#)

Parameter	Details
Product Name	POGLUT1 Rabbit pAb
Remarks/Alias	Rumi; CLP46; MDSRP; C3orf9; KTELC1; LGMD2Z; MDS010; hCLP46; KDELC1; LGMDR21; POGLUT1
Species	Human
Gene ID	56983
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 160-280 of human POGLUT1 (NP_689518.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8NBL1
Protein Weight	46kDa
Shipping	Ice bag