


Anti-Pold1 Monoclonal Antibody For Wb, Elisa - P28340

Item Number: CM0009732



Introduction

Rabbit monoclonal antibody targeting human DNA polymerase delta catalytic subunit (POLD1), validated for Western blot (WB) and ELISA. Cross-reacts with human, mouse, and rat samples. Ideal for DNA replication and repair studies.

[Learn More](#)

Parameter	Value
Product Name	POLD1 Rabbit mAb
Remarks/Alias	CDC2; MDPL; POLD; CRC510; POLD1
Species	Human
GeneID (Human)	5424
GeneID	5424
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1000-1107 of human POLD1 (P28340).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28340
Protein Weight	124kDa
Shipping	Ice bag