


Anti-Pold2 Polyclonal Antibody For Wb - P49005

Item Number: CM0030808



Introduction

Anti-POLD2 rabbit polyclonal antibody validated for WB and ELISA, detecting human POLD2 (DNA polymerase delta subunit 2) with cross-reactivity to mouse and rat; widely used in DNA replication and repair research.

[Learn More](#)

Parameter	Value
Product Name	POLD2 Rabbit pAb
Remarks/Alias	POLD2
Species	Human
GeneID (Human)	5425
GeneID	5425
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 120-360 of human POLD2 (NP_006221.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49005
Protein Weight	51kDa
Shipping	Ice bag