


Anti-Pold3 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15054

Item Number: CM0026845



Introduction

POLD3 Rabbit pAb, a rabbit polyclonal antibody targeting DNA polymerase delta subunit 3 for WB, IHC-P, IF/ICC, ELISA; suitable for human and mouse samples. Research use in DNA replication and repair.

[Learn More](#)

Parameter	Value
Product Name	POLD3 Rabbit pAb
Remarks/Alias	P66; P68; PPP1R128; POLD3
Species	Human
GeneID (Human)	10714
GeneID	10714
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 197-466 of human POLD3 (NP_006582.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15054
Protein Weight	51kDa
Shipping	Ice bag