

Anti-Poldip2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y2S7

Item Number: CM0015476



Introduction

This polyclonal antibody is raised in rabbit against human POLDIP2 and validated for WB and ELISA. It recognizes the 42 kDa protein and cross-reacts with mouse and rat orthologs, making it a valuable tool for DNA repair studies.

[Learn More](#)

Parameter	Value
Product Name	POLDIP2 Rabbit pAb
Remarks/Alias	p38; POLD4; PDIP38; POLDIP2
Species	Human
Gene ID (Human)	26073
Gene ID	26073
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 250-350 of human POLDIP2 (NP_056399.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2S7
Protein Weight	42kDa
Shipping	Ice bag