

Anti-Pole2 Rabbit Polyclonal Antibody For Wb And Elisa - P56282

Item Number: CM0017443



Introduction

Rabbit anti-POLE2 polyclonal antibody validated for Western blot and ELISA; reacts with human and mouse. DNA polymerase epsilon subunit 2 involved in DNA replication and repair.

[Learn More](#)

Parameter	Value
Product Name	POLE2 Rabbit pAb
Remarks/Alias	DPE2; POLE2
Species	Human
GeneID (Human)	5427
GeneID	5427
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P56282
Protein Weight	60kDa
Shipping	Ice bag