

Anti-Pole4 Polyclonal Antibody For Wb, Elisa - Q9Nr33

Item Number: CM0026901

Introduction

Rabbit polyclonal antibody targeting human POLE4 (p12), a DNA polymerase epsilon accessory subunit. Validated in WB and ELISA, ideal for DNA replication and repair research.

[Learn More](#)

Parameter	Value
Product Name	POLE4 Rabbit pAb
Remarks/Alias	p12; YHHQ1; POLE4
Species	Human
Gene ID	56655
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-117 of human POLE4 (NP_063949.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NR33
Protein Weight	12kDa
Shipping	Ice bag