

Anti-Pole4 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nr33

Item Number: CM0023859

Introduction

Rabbit polyclonal antibody targeting human POLE4, a DNA polymerase epsilon subunit involved in DNA replication and repair. Validated for Western blot and ELISA; cross-reacts with human and mouse. Ideal for chromatin studies.

[Learn More](#)

Parameter	Value
Product Name	POLE4 Rabbit pAb
Remarks/Alias	POLE4; YHHQ1; p12; DNA polymerase epsilon subunit 4
Species	Human
Gene ID	56655
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human POLE4 (NP_063949.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NR33
Protein Weight	12kDa
Shipping	Ice bag