

Anti-PolG2 Rabbit Polyclonal Antibody - Q9Uhn1

Item Number: CM0027994



Introduction

High-quality anti-POLG2 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Specific to human, mouse, rat POLG2 (Q9UHN1), mitochondrial DNA polymerase accessory subunit.

[Learn More](#)

Parameter	Value
Product Name	POLG2 Rabbit pAb
Remarks/Alias	HP55; POLB; PEOA4; POLGB; MTPOLB; MTDPS16; MTDPS16A; MTDPS16B; POLG-BETA; POLG2
Species	Human
GeneID (human)	11232
GeneID	11232
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 216-485 of human POLG2 (NP_009146.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHN1
Protein Weight	55kDa
Shipping	Ice bag