


Anti-POLL Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ugp5

Item Number: CM0030455



Introduction

Rabbit polyclonal antibody against human DNA polymerase lambda (POLL), validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for research in DNA repair pathways including base excision repair and non-homologous end joining.

[Learn More](#)

Parameter	Value
Product Name	POLL Rabbit pAb
Remarks/Alias	BETAN; POLKAPPA; POLL
Species	Human
GeneID (Human)	27343
GeneID	27343
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 421-575 of human POLL (NP_037406.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UGP5
Protein Weight	63kDa
Shipping	Ice bag