


Anti-Polr1B Polyclonal Antibody For Wb, Elisa - Q9H9Y6

Item Number: CM0018514



Introduction

Rabbit polyclonal antibody targeting POLR1B (RPA135), a core subunit of RNA polymerase I. Validated for WB and ELISA, cross-reacts with Human and Mouse. Ideal for nucleolar transcription studies.

[Learn More](#)

Parameter	Value
Product Name	POLR1B Rabbit pAb
Remarks/Alias	A135; RPA2; TCS4; RPA135; Rpo1-2; POLR1B
Species	Human
Gene ID (Human)	84172
Gene ID	84172
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 670-900 of human POLR1B (NP_061887.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H9Y6
Protein Weight	128kDa
Shipping	Ice bag