


## Anti-Polr1C Rabbit Polyclonal Antibody For Wb And Elisa - O15160

Item Number: CM0014088



### Introduction

Rabbit polyclonal antibody targeting human POLR1C (RPAC1), validated for Western blot and ELISA. Essential for studying RNA polymerase I/III transcription, nucleolar function, and leukodystrophy-related pathways. - O15160

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | POLR1C Rabbit pAb                                                                                                                       |
| Remarks/Alias    | AC40; RPA5; TCS3; HLD11; RPA39; RPA40; RPAC1; RPC40; POLR1C                                                                             |
| Species          | Human                                                                                                                                   |
| GeneID (Human)   | 9533                                                                                                                                    |
| GeneID           | 9533                                                                                                                                    |
| Immunogen        | Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human POLR1C (NP_976035.1). |
| Source           | Rabbit                                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                                   |
| Application      | WB, ELISA                                                                                                                               |
| Cross Reactivity | Human                                                                                                                                   |
| SWISS            | O15160                                                                                                                                  |
| Protein Weight   | 39kDa                                                                                                                                   |
| Shipping         | Ice bag                                                                                                                                 |