


Anti-Polr1E Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q9Gzs1

Item Number: CM0013738



Introduction

High-quality rabbit polyclonal antibody against POLR1E (PAF53/RPA49), a Pol I subunit vital for rRNA transcription. Validated in WB, IP, ELISA; reacts with human, mouse, rat. Ideal for nucleolar and transcription research.

[Learn More](#)

Parameter	Value
Product Name	POLR1E Rabbit pAb
Remarks/Alias	A49; PAF53; PRAF1; RPA49; POLR1E
Species	Human
GeneID (Human)	64425
GeneID	64425
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-419 of human POLR1E (NP_071935.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZS1
Protein Weight	47kDa
Shipping	Ice bag