


Anti-Polr2A (N-Terminal) Rabbit Mab For Wb, Ip, Chip, Elisa - P24928

Item Number: CM0002675



Introduction

High-quality anti-POLR2A (N-terminal) rabbit monoclonal antibody validated for WB, IP, CHIP, and ELISA; cross-reacts with human, mouse, and rat. Suitable for RNA polymerase II transcription research.

[Learn More](#)

Parameter	Value
Product Name	POLR2A (N-terminal) Rabbit mAb
Remarks/Alias	RPB1; RPO2; POLR2; POLRA; RPBh1; RPOL2; NEDHIB; RplILS; hsRPB1; hRPB220; POLR2A (N-terminal)
Species	Human
GeneID (Human)	5430
GeneID	5430
Immunogen	Synthetic peptide corresponding to amino acids 550-650 of human POLR2A NP_000928.1
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, CHIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24928
Protein Weight	217 kDa
Shipping	Ice bag