


Anti-Polr2B Polyclonal Antibody For Wb, Ihc-P, Elisa - P30876

Item Number: CM0030289



Introduction

Anti-POLR2B rabbit polyclonal antibody validated for WB, IHC-P, ELISA applications. Recognizes human, mouse, rat POLR2B (134kDa). Recombinant immunogen (aa 1-260). Essential for RNA polymerase II transcription and chromatin research.

[Learn More](#)

Parameter	Specification
Product Name	POLR2B Rabbit pAb
Remarks/Alias	RPB2, POL2RB, hRPB140, POLR2B
Species	Human
GeneID (human)	5431
GeneID	5431
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human POLR2B (NP_000929.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30876
Protein Weight	134kDa
Shipping	Ice bag