

Anti-Polr2D Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - O15514

Item Number: CM0020060



Introduction

Rabbit polyclonal antibody targeting human POLR2D, a core RNA polymerase II subunit. Validated for WB, IF/ICC, IP, ELISA; cross-reacts with mouse. Ideal for transcription research.

[Learn More](#)

Parameter	Value
Product Name	POLR2D Rabbit pAb
Remarks/Alias	RBP4; RPB4; RPB16; HSRBP4; HSRPB4; POLR2D
Species	Human
Gene ID	5433
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-142 of human POLR2D (NP_004796.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15514
Protein Weight	16kDa
Shipping	Ice bag